Low Window[®]

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/ Bibliografia

- **Zaniol T, Zaniol A.**
La tecnica Low Window nel grande rialzo di seno mascellare.
Implant Tribune 2018 March.
- **Zaniol T, Zaniol A.**
A Rational Approach to Sinus Augmentation: The Low Window Sinus Lift.
Case Rep Dent. 2017;2017:7610607.
- **Arosio P, Greco GB, Zaniol T, Iezzi G, Perrotti V, Di Stefano DA.**
Sinus augmentation and concomitant implant placement in low bone-density sites.
A retrospective study on an undersized drilling protocol and primary stability.
Clin Implant Dent Relat Res. 2017 Nov 22. [Epub ahead of print].
- **Mandelaris GA, Rosenfeld AL.**
A novel approach to the antral sinus bone graft technique: the use of a prototype cutting guide for precise outlining of the lateral wall. A case report.
Int J Periodontics Restorative Dent 2008;28:569-575.
- **Mandelaris GA, Rosenfeld AL.** *Alternative applications of guided surgery: precise outlining of the lateral window in antral sinus bone grafting.*
J Oral Maxillofac Surg 2009;67:23-30.
- **Stern A, Green J.**
Sinus lift procedures: an overview of current techniques.
Dent Clin North Am. 2012 Jan;56(1):219-33.
- **Guerrero JS, Al-Jandan BA.**
Lateral wall sinus floor elevation for implant placement: revisiting fundamentals and the surgical technique.
J Calif Dent Assoc 2013;41:185-187,190-195.
- **Cho SC, Wallace SS, Froum SJ, Tarnow DP.**
Influence of anatomy on Schneiderian membrane perforations during sinus elevation surgery: three-dimensional analysis.
Pract Proced Aesthet Dent 2001;13:160-163.
- **Vercellotti T, De Paoli S, Nevins M.**
The piezoelectric bony window osteotomy and sinus membrane elevation: introduction of a new technique for simplification of the sinus augmentation procedure.
Int J Periodontics Restorative Dent. 2001 Dec;21(6):561-7.
- **Velasco-Torres M, Padial-Molina M, Alarcón JA, O'Valle F, Catena A, Galindo-Moreno P.**
Maxillary Sinus Dimensions With Respect to the Posterior Superior Alveolar Artery Decrease With Tooth Loss.
Implant Dent. 2016 Aug;25(4):464-70.
- **Di Stefano DA, Cazzaniga A.**
Tecniche rigenerative e ricostruttive in chirurgia pre- e perimplantare.
Elsevier. Milano, 2012.